

|   |   |  |           |
|---|---|--|-----------|
| <b>IHE Integration Statement</b>  |   | <b>Date:</b>   | May, 2006 |
| <b>Vendor</b>   | <b>Product Name</b>   | <b>Version</b>   |           |
| Eastman Kodak Company   | <i>Kodak Versatile Intelligent Patient Archive / VIParchive</i>   | 10.0   |           |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below: |   |  |           |
| <b>Integration Profiles Implemented</b>   | <b>Actors Implemented</b>   | <b>Options Implemented</b>   |           |
| Cross-Enterprise Document Sharing (XDS)   | Document Repository<br><br>Document Source, grouped with the Document Repository<br><br>Note: The "Register Document" transaction is sent directly, without the "Provide and Register Document"   | None<br><br>Multiple Document Submission<br>Document Life Cycle Management<br>Folder Management  |           |
| Retrieve Information for Display (RID)  | Information Source<br><br>Notes:<br><ul style="list-style-type: none"> <li>• GMT times are not supported in the lowerDateTime and upperDateTime parameters</li> <li>• Format conversions are not supported</li> <li>• DICOM Studies are also made accessible through this transaction, in zip or multipart/related format.</li> <li>• Private extensions are available, as a part of the VIParchive API.</li> </ul> | Persistent Document<br>Summary of All Reports<br>Summary of Laboratory Reports<br>Summary of Radiology Reports<br>Summary of Cardiology Reports<br>Summary of Surgery Reports<br>Summary of Intensive Care Reports<br>Summary of Emergency Reports<br>Summary of Discharge Reports<br>Summary of Prescriptions |           |
| Audit Trail and Node Authentication (ATNA)  | Secure Node (for the XDS and RID actors)<br>Document Repository<br>Embedded Repository<br>Information Source  | None   |           |
| Consistent Time (CT)  | Time Client   | None   |           |
| Teaching file and Clinical trial Export (TCE)   | Receiver  | Additional Teaching File Information   |           |
| <b>Internet address for Kodak IHE information:</b> <a href="http://www.kodak.com/go/ihe">www.kodak.com/go/ihe</a>   |   |  |           |
| <b>Links to Standards Conformance Statements for the Implementation</b>   |   |  |           |
| <b>DICOM</b>  | <a href="http://www.kodak.com/go/dicom">www.kodak.com/go/dicom</a>  |  |           |
| <b>Links to general information on IHE</b>  |   |  |           |
| North America:<br><a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>  | Europe:<br><a href="http://www.ihe-europe.org">www.ihe-europe.org</a>   | Japan:<br><a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>   |           |